



KIT CONTAINS:

- N.2 VACUUM TUBES WITH SEPARATOR GEL AND ANTICOAGULANT

ACCESSORIES SET CONTAINS:

- BUTTERFLY NEEDLE, TOURNIQUET
- PACTH, VACUAINER
- STERILE LUER-LOCK SYRINGE 5ml WITH NEEDLE
- STERILE LUER-LOCK SYRINGE 1ml WITH NEEDLE
- N.1 INSTRUCTION USE AND MSDN

NB.:

1) CENTRIFUGE WITH FIXED ROTOR:
For the centrifuge recommended by the manufacturer, the number of revolutions is rpm 4000 rpm for 10 minutes.

2) CENTRIFUGE WITH ADJUSTABLE ROTOR:
Please contact the manufacturer of the kit to calculate the conversion from rcf to rpm or using the following conversion formula:

$$RCF = \left(\frac{RPM}{1000} \right)^2 \times r \times 1,118 \Rightarrow RPM = \sqrt{\frac{RCF}{r \times 1,118}} \times 1000$$

REF GPH1108D

8 043400 003969

